

Tissue alpha-L-fucosidase

Catalog No: #AP83336



Package Size: #AP83336-1 50ug #AP83336-2 100ug #AP83336-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tissue alpha-L-fucosidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Alpha-L-fucosidase I,Alpha-L-fucoside fucohydrolase 1,Alpha-L-fucosidase 1,Fuca
Accession No.	Uniprot ID:Q99LJ1
Uniprot	Q99LJ1
GeneID	71665;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fuca1

Note: This product is for in vitro research use only