

L-gulonolactone oxidase

Catalog No: #AP83338



Package Size: #AP83338-1 50ug #AP83338-2 100ug #AP83338-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	L-gulonolactone oxidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	LGO,L-gulono-gamma-lactone oxidase,GLO
Accession No.	Uniprot ID:Q8HXW0
Uniprot	Q8HXW0
GeneID	396759;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GULO

Note: This product is for in vitro research use only