

Gamma-interferon-inducible lysosomal thiol reductase

Catalog No: #AP83339



Package Size: #AP83339-1 50ug #AP83339-2 100ug #AP83339-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Gamma-interferon-inducible lysosomal thiol reductase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Gamma-interferon-inducible protein IP-30,Lysosomal thiol reductase IP30,Gilt,Ip30
Accession No.	Uniprot ID:Q9ESY9
Uniprot	Q9ESY9
GeneID	65972;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ifi30

Note: This product is for in vitro research use only