

GDNF family receptor alpha-1

Catalog No: #AP83344



Package Size: #AP83344-1 50ug #AP83344-2 100ug #AP83344-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GDNF family receptor alpha-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	GDNF receptor alpha-1, TGF-beta-related neurotrophic factor receptor 1, GDNFRA
Accession No.	Uniprot ID:O13156
Uniprot	O13156
GeneID	395994;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: GFRA1

Note: This product is for in vitro research use only