

Protein fem-1 homolog A

Catalog No: #AP83346



Package Size: #AP83346-1 50ug #AP83346-2 100ug #AP83346-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein fem-1 homolog A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	FEM1a,FEM1-alpha
Accession No.	Uniprot ID:Q29RM5
Uniprot	Q29RM5
GeneID	516550;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FEM1A

Note: This product is for in vitro research use only