

Ferritin light chain 1

Catalog No: #AP83351



Package Size: #AP83351-1 50ug #AP83351-2 100ug #AP83351-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ferritin light chain 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Ferritin L subunit 1,Ftl
Accession No.	Uniprot ID:P02793
Uniprot	P02793
GeneID	100360087;29292;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ftl1

Note: This product is for in vitro research use only