

Follistatin-related protein 3

Catalog No: #AP83352



Package Size: #AP83352-1 50ug #AP83352-2 100ug #AP83352-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Follistatin-related protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Follistatin-like protein 3
Accession No.	Uniprot ID:Q1LZB9
Uniprot	Q1LZB9
GeneID	515367;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: FSTL3

Note: This product is for in vitro research use only