

Protein FAM3B

Catalog No: #AP83355



Package Size: #AP83355-1 50ug #AP83355-2 100ug #AP83355-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein FAM3B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cytokine-like protein 2-21,Pancreatic-derived factor,PANDER,C21orf11,C21orf76,PRED44,UNQ320/PRO365
Accession No.	Uniprot ID:P58499
Uniprot	P58499
GeneID	54097;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FAM3B

Note: This product is for in vitro research use only