

Growth hormone receptor

Catalog No: #AP83360



Package Size: #AP83360-1 50ug #AP83360-2 100ug #AP83360-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Growth hormone receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	GH receptor,Somatotropin receptor
Accession No.	Uniprot ID:Q02092
Uniprot	Q02092
GeneID	408184;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GHR

Note: This product is for in vitro research use only