

Embryonic stem cell-specific 5-hydroxymethylcytosine-binding protein

Catalog No: #AP83364

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP83364-1 50ug #AP83364-2 100ug #AP83364-3 1mg

Description

Product Name	Embryonic stem cell-specific 5-hydroxymethylcytosine-binding protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ES cell-specific 5hmC-binding protein,Putative peptidase SRAPD1,SRAP domain-containing protein 1,Srapd1
Accession No.	Uniprot ID:Q8R1M0
Uniprot	Q8R1M0
GeneID	232210;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hmces

Note: This product is for in vitro research use only