

Tissue alpha-L-fucosidase

Catalog No: #AP83367



Package Size: #AP83367-1 50ug #AP83367-2 100ug #AP83367-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tissue alpha-L-fucosidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Alpha-L-fucosidase I,Alpha-L-fucoside fucohydrolase 1,Alpha-L-fucosidase 1
Accession No.	Uniprot ID:Q2KIM0
Uniprot	Q2KIM0
GeneID	509522;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FUCA1

Note: This product is for in vitro research use only