

Follistatin-related protein 1

Catalog No: #AP83370

Package Size: #AP83370-1 50ug #AP83370-2 100ug #AP83370-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Follistatin-related protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Follistatin-like protein 1
Accession No.	Uniprot ID:Q58D84
Uniprot	Q58D84
GeneID	534482;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: FSTL1

Note: This product is for in vitro research use only