

Endoplasmin

Catalog No: #AP83371



Package Size: #AP83371-1 50ug #AP83371-2 100ug #AP83371-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Endoplasmin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Heat shock 108 kDa protein,HSP 108,Heat shock protein 90 kDa beta member 1,Transferrin-binding protein,TRA1
Accession No.	Uniprot ID:P08110
Uniprot	P08110
GeneID	374163;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HSP90B1

Note: This product is for in vitro research use only