

EGF, latrophilin and seven transmembrane domain-containing protein 1

Catalog No: #AP83375

Package Size: #AP83375-1 50ug #AP83375-2 100ug #AP83375-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EGF, latrophilin and seven transmembrane domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	EGF-TM7-latrophilin-related protein,ETL protein,Etl
Accession No.	Uniprot ID:Q9ESC1
Uniprot	Q9ESC1
GeneID	64124;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Eltd1

Note: This product is for in vitro research use only