

Glutamate receptor 2

Catalog No: #AP83380



Package Size: #AP83380-1 50ug #AP83380-2 100ug #AP83380-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glutamate receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GluR-2,AMPA-selective glutamate receptor 2,GluR-B,GluR-K2,Glutamate receptor ionotropic, AMPA 2,GluA2,GLUR2
Accession No.	Uniprot ID:P42262
Uniprot	P42262
GeneID	2891;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GRIA2

Note: This product is for in vitro research use only