

## Gasdermin-C

Catalog No: #AP83381



Package Size: #AP83381-1 50ug #AP83381-2 100ug #AP83381-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Gasdermin-C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Gasdermin-C1, Melanoma-derived leucine zipper-containing extranuclear factor, mMLZE, Gsdmc1, Mlze
Accession No.	Uniprot ID:Q99NB5
Uniprot	Q99NB5
GeneID	83492;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Gsdmc

Note: This product is for in vitro research use only