

Gamma-aminobutyric acid type B receptor subunit 1

Catalog No: #AP83393



Package Size: #AP83393-1 50ug #AP83393-2 100ug #AP83393-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Gamma-aminobutyric acid type B receptor subunit 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	GABA-B receptor 1
Accession No.	Uniprot ID:Q9Z0U4
Uniprot	Q9Z0U4
GeneID	81657;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Gabbr1

Note: This product is for in vitro research use only