

High mobility group protein HMGI-C

Catalog No: #AP83395



Package Size: #AP83395-1 50ug #AP83395-2 100ug #AP83395-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	High mobility group protein HMGI-C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	High mobility group AT-hook protein 2,HMGIC
Accession No.	Uniprot ID:P52926
Uniprot	P52926
GeneID	8091;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HMGA2

Note: This product is for in vitro research use only