

## Filamin-A

Catalog No: #AP83400



Package Size: #AP83400-1 50ug #AP83400-2 100ug #AP83400-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Filamin-A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FLN-A,Actin-binding protein 280,ABP-280,Alpha-filamin,Endothelial actin-binding protein,Filamin-1,Non-muscle filamin,FLN,FLN1
Accession No.	Uniprot ID:P21333
Uniprot	P21333
GeneID	2316;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:FLNA

Note: This product is for in vitro research use only