

Platelet glycoprotein Ib alpha chain

Catalog No: #AP83405



Package Size: #AP83405-1 50ug #AP83405-2 100ug #AP83405-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Platelet glycoprotein Ib alpha chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GP-Ib alpha, Antigen CD42b-alpha, CD42b
Accession No.	Uniprot ID:P07359
Uniprot	P07359
GeneID	2811;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: GP1BA

Note: This product is for in vitro research use only