

Follistatin-related protein 1

Catalog No: #AP83406



Package Size: #AP83406-1 50ug #AP83406-2 100ug #AP83406-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Follistatin-related protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Follistatin-like protein 1,FRP
Accession No.	Uniprot ID:Q12841
Uniprot	Q12841
GeneID	11167;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FSTL1

Note: This product is for in vitro research use only