

Glycogen synthase kinase-3 alpha

Catalog No: #AP83411



Package Size: #AP83411-1 50ug #AP83411-2 100ug #AP83411-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glycogen synthase kinase-3 alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	GSK-3 alpha, Serine/threonine-protein kinase GSK3A
Accession No.	Uniprot ID:Q2NL51
Uniprot	Q2NL51
GeneID	606496;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Gsk3a

Note: This product is for in vitro research use only