

Coagulation factor VII

Catalog No: #AP83413



Package Size: #AP83413-1 50ug #AP83413-2 100ug #AP83413-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Coagulation factor VII
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Serum prothrombin conversion accelerator
Accession No.	Uniprot ID:Q8K3U6
Uniprot	Q8K3U6
GeneID	260320;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: F7

Note: This product is for in vitro research use only