

Gastrotropin

Catalog No: #AP83415



Package Size: #AP83415-1 50ug #AP83415-2 100ug #AP83415-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Gastrotropin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GT,Fatty acid-binding protein 6,Ileal lipid-binding protein,ILBP,Intestinal 15 kDa protein,I-15P,Intestinal bile acid-binding protein,I-BABP,ILBP,ILLBP
Accession No.	Uniprot ID:P51161
Uniprot	P51161
GeneID	2172;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FABP6

Note: This product is for in vitro research use only