

Flotillin-2

Catalog No: #AP83416



Package Size: #AP83416-1 50ug #AP83416-2 100ug #AP83416-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Flotillin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Epidermal surface antigen,ESA,Membrane component chromosome 17 surface marker 1,ESA1,M17S1
Accession No.	Uniprot ID:Q14254
Uniprot	Q14254
GeneID	2319;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FLOT2

Note: This product is for in vitro research use only