

Forkhead box protein O1

Catalog No: #AP83417



Package Size: #AP83417-1 50ug #AP83417-2 100ug #AP83417-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Forkhead box protein O1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Forkhead box protein O1A,Forkhead in rhabdomyosarcoma,Fkhr,Foxo1a
Accession No.	Uniprot ID:Q9R1E0
Uniprot	Q9R1E0
GeneID	56458;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Foxo1

Note: This product is for in vitro research use only