

# Follitropin subunit beta

Catalog No: #AP83418



Package Size: #AP83418-1 50ug #AP83418-2 100ug #AP83418-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Follitropin subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Follicle-stimulating hormone beta subunit,FSH-B,Follitropin beta chain
Accession No.	Uniprot ID:P01228
Uniprot	P01228
GeneID	396895;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:FSHB

Note: This product is for in vitro research use only