

## PASK (Phospho-Thr1165) Antibody

Catalog No: #12526

Package Size: #12526-1 50ul #12526-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	PASK (Phospho-Thr1165) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	PASK (Phospho-Thr1165) Antibody detects endogenous levels of PASK only when phosphorylated at Thr1165
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human PASK (Phospho-Thr1165)
Target Name	PASK
Modification	Phospho
Other Names	PASK, HPASK, KIAA0135, PAS Kinase, PAS-kinase, STK37, PASKIN
Accession No.	Swiss-Prot#: Q96RG2NCBI Gene ID: 23178
Uniprot	Q96RG2
GeneID	23178;
Target Species	human
Calculated MW	142kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000

## Images



Western Blot analysis of lysates of 1. Jurkat cell and 2. Jurkat cell treated with 100ng/mL LPS for 30min, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only