

Gastrin/cholecystokinin type B receptor

Catalog No: #AP83421



Package Size: #AP83421-1 50ug #AP83421-2 100ug #AP83421-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Gastrin/cholecystokinin type B receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CCK-B receptor,Cholecystokinin-2 receptor,CCK2-R
Accession No.	Uniprot ID:P56481
Uniprot	P56481
GeneID	12426;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cckbr

Note: This product is for in vitro research use only