

Guanylate-binding protein 4

Catalog No: #AP83424



Package Size: #AP83424-1 50ug #AP83424-2 100ug #AP83424-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Guanylate-binding protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GTP-binding protein 4,GBP-4,Guanine nucleotide-binding protein 4
Accession No.	Uniprot ID:Q96PP9
Uniprot	Q96PP9
GeneID	115361;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GBP4

Note: This product is for in vitro research use only