

L-gulonolactone oxidase

Catalog No: #AP83429



Package Size: #AP83429-1 50ug #AP83429-2 100ug #AP83429-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	L-gulonolactone oxidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	LGO,L-gulono-gamma-lactone oxidase,GLO
Accession No.	Uniprot ID:P58710
Uniprot	P58710
GeneID	268756;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gulo

Note: This product is for in vitro research use only