

F-box only protein 32

Catalog No: #AP83430



Package Size: #AP83430-1 50ug #AP83430-2 100ug #AP83430-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--------------------------------------------------------|
| Product Name | F-box only protein 32 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | Atrogin-1, Muscle atrophy F-box protein, MAFbx |
| Accession No. | Uniprot ID:Q91Z62 |
| Uniprot | Q91Z62 |
| GeneID | 171043; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Fbxo32

Note: This product is for in vitro research use only