

IgG receptor FcRn large subunit p51

Catalog No: #AP83431



Package Size: #AP83431-1 50ug #AP83431-2 100ug #AP83431-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	IgG receptor FcRn large subunit p51
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	FcRn,IgG Fc fragment receptor transporter alpha chain,Neonatal Fc receptor,FcRn
Accession No.	Uniprot ID:Q61559
Uniprot	Q61559
GeneID	14132;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fcgrt

Note: This product is for in vitro research use only