

Free fatty acid receptor 4

Catalog No: #AP83432



Package Size: #AP83432-1 50ug #AP83432-2 100ug #AP83432-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Free fatty acid receptor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	G-protein coupled receptor 120,G-protein coupled receptor 129,G-protein coupled receptor GT01,G-protein coupled receptor PGR4,Omega-3 fatty acid receptor 1,GPR120,GPR129,O3FAR1,PGR4
Accession No.	Uniprot ID:Q5NUL3
Uniprot	Q5NUL3
GeneID	338557;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FFAR4

Note: This product is for in vitro research use only