

Fibroblast growth factor 3

Catalog No: #AP83434



Package Size: #AP83434-1 50ug #AP83434-2 100ug #AP83434-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fibroblast growth factor 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	FGF-3,Heparin-binding growth factor 3,HBGF-3,Proto-oncogene Int-2,Fgf-3,Int-2
Accession No.	Uniprot ID:P05524
Uniprot	P05524
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fgf3

Note: This product is for in vitro research use only