

## Receptor-type tyrosine-protein kinase FLT3

Catalog No: #AP83435



Package Size: #AP83435-1 50ug #AP83435-2 100ug #AP83435-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Receptor-type tyrosine-protein kinase FLT3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	FL cytokine receptor, Fetal liver kinase 2, FLK-2, Fms-like tyrosine kinase 3, FLT-3, Tyrosine-protein kinase receptor flk-2, CD135, Flk-2, Flt-3
Accession No.	Uniprot ID:Q00342
Uniprot	Q00342
GeneID	14255;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Flt3

Note: This product is for in vitro research use only