

Fibrinogen C domain-containing protein 1

Catalog No: #AP83437



Package Size: #AP83437-1 50ug #AP83437-2 100ug #AP83437-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fibrinogen C domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	UNQ701/PRO1346
Accession No.	Uniprot ID:Q8N539
Uniprot	Q8N539
GeneID	84929;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: FIBCD1

Note: This product is for in vitro research use only