

Protein fantom

Catalog No: #AP83439



Package Size: #AP83439-1 50ug #AP83439-2 100ug #AP83439-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein fantom
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Nephrocystin-8,RPGR-interacting protein 1-like protein,RPGRIP1-like protein,FTM,KIAA1005,NPHP8
Accession No.	Uniprot ID:Q68CZ1
Uniprot	Q68CZ1
GeneID	23322;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RPGRIP1L

Note: This product is for in vitro research use only