

## Follistatin-related protein 3

Catalog No: #AP83444



Package Size: #AP83444-1 50ug #AP83444-2 100ug #AP83444-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Follistatin-related protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Follistatin-like protein 3,Follistatin-related gene protein,FLRG,UNQ674/PRO1308
Accession No.	Uniprot ID:O95633
Uniprot	O95633
GeneID	10272;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:FSTL3

Note: This product is for in vitro research use only