

Guanine nucleotide-binding protein G(s) subunit alpha

Catalog No: #AP83447



Package Size: #AP83447-1 50ug #AP83447-2 100ug #AP83447-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Guanine nucleotide-binding protein G(s) subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Adenylate cyclase-stimulating G alpha protein,GNAS1
Accession No.	Uniprot ID:P63091
Uniprot	P63091
GeneID	403943;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GNAS

Note: This product is for in vitro research use only