

## Gypsy retrotransposon integrase-like protein 1

Catalog No: #AP83450



Package Size: #AP83450-1 50ug #AP83450-2 100ug #AP83450-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Gypsy retrotransposon integrase-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GIN-1,Ty3/Gypsy integrase 1,Zinc finger H2C2 domain-containing protein,TGIN1,ZH2C2
Accession No.	Uniprot ID:Q9NXP7
Uniprot	Q9NXP7
GeneID	54826;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:GIN1

Note: This product is for in vitro research use only