

Gonadotropin-releasing hormone receptor

Catalog No: #AP83454



Package Size: #AP83454-1 50ug #AP83454-2 100ug #AP83454-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Gonadotropin-releasing hormone receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	GnRH receptor, Luteinizing hormone-releasing hormone receptor, LHRH
Accession No.	Uniprot ID:P49922
Uniprot	P49922
GeneID	397515;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: GNRHR

Note: This product is for in vitro research use only