

Galactocerebrosidase

Catalog No: #AP83456



Package Size: #AP83456-1 50ug #AP83456-2 100ug #AP83456-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Galactocerebrosidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	GALCERase, Galactocerebroside beta-galactosidase, Galactosylceramidase, Galactosylceramide beta-galactosidase
Accession No.	Uniprot ID:P54818
Uniprot	P54818
GeneID	14420;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Galc

Note: This product is for in vitro research use only