

Gamma-aminobutyric acid receptor subunit alpha-3

Catalog No: #AP83458



Package Size: #AP83458-1 50ug #AP83458-2 100ug #AP83458-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Gamma-aminobutyric acid receptor subunit alpha-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	GABA(A) receptor subunit alpha-3,Gabra-3
Accession No.	Uniprot ID:P26049
Uniprot	P26049
GeneID	14396;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gabra3

Note: This product is for in vitro research use only