

Guanine nucleotide-binding protein G(s) subunit alpha isoforms short



Catalog No: #AP83459

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP83459-1 50ug #AP83459-2 100ug #AP83459-3 1mg

Description

Product Name	Guanine nucleotide-binding protein G(s) subunit alpha isoforms short
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Adenylate cyclase-stimulating G alpha protein,GNAS1
Accession No.	Uniprot ID:P04896
Uniprot	O18979
GeneID	281793;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GNAS

Note: This product is for in vitro research use only