

## FAS-associated factor 2

Catalog No: #AP83461



Package Size: #AP83461-1 50ug #AP83461-2 100ug #AP83461-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	FAS-associated factor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Protein ETEA,UBX domain-containing protein 3B,UBX domain-containing protein 8,ETEA,KIAA0887,UBXD8,UBXN3B
Accession No.	Uniprot ID:Q96CS3
Uniprot	Q96CS3
GeneID	23197;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:FAF2

Note: This product is for in vitro research use only