

Frizzled-4

Catalog No: #AP83463



Package Size: #AP83463-1 50ug #AP83463-2 100ug #AP83463-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Frizzled-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Fz-4, FZ4
Accession No.	Uniprot ID:Q9IA05
Uniprot	Q9IA05
GeneID	373887;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: FZD4

Note: This product is for in vitro research use only