

N-acetylgalactosamine-6-sulfatase

Catalog No: #AP83465



Package Size: #AP83465-1 50ug #AP83465-2 100ug #AP83465-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	N-acetylgalactosamine-6-sulfatase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Chondroitinsulfatase,Chondroitinase,Galactose-6-sulfate sulfatase,N-acetylgalactosamine-6-sulfate sulfatase,GalNAc6S sulfatase
Accession No.	Uniprot ID:Q32KJ6
Uniprot	Q32KJ6
GeneID	292073;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Galns

Note: This product is for in vitro research use only