

## Glucagon-like peptide 1 receptor

Catalog No: #AP83468



Package Size: #AP83468-1 50ug #AP83468-2 100ug #AP83468-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Glucagon-like peptide 1 receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	GLP-1 receptor, Glpr
Accession No.	Uniprot ID:P32301
Uniprot	P32301
GeneID	25051;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Glp1r

Note: This product is for in vitro research use only