

Follitropin subunit beta

Catalog No: #AP83470



Package Size: #AP83470-1 50ug #AP83470-2 100ug #AP83470-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Follitropin subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Follicle-stimulating hormone beta subunit,FSH-B,Follitropin beta chain
Accession No.	Uniprot ID:Q60687
Uniprot	Q60687
GeneID	14308;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fshb

Note: This product is for in vitro research use only